
YeoMed

Smart Baby Scale



Model: Y-11PC







Functions



Body Length



Weight



WHO Child Growth Standard




Head Circumference



Chest Circumference



Head-to-Rump Length



- 10.1" Touch Display
- 2 Magnetic Card Reader
- 3 Bluetooth ruler
- 4 Code Scanner
- Double Ultrasonic Height Sensor

- 6 Length Ruler
- USB*2; RJ45 for data transmit
- Movable skateboard for Length measurement
- Weighing Pan
- Foot Pads

Operation Instructions





Start Measurement



Cloud Server



Measurement Completed



Weight, Body Length



Head Circumference, Chest Circumference, Head-to-Rump Length, WHO Child Growth Standard

Features


























Data Interconnection



Technical Parameters

$T \sim \sim$	hnion	l Parameters
IEC	IIIIGa	i Farailleters

Model	Y-11PC	
Login Methods	QR Code/Barcode; Magnetic Card; Phone Number	
Login Devices	Code Scanner	QR Code, Barcode
	Magnetic Card Reader	Supports ISO/IBM/CA/DMV/AAMVA formats
	Processor	Four-core processor
	Memory	2G
	Display	10.1 inch 800*1280 multi-touch
	Storage	32G
	Operating System	Android 11
Device	Speaker	8Ω, 2W
Parameters	Height Measurement Range	30-105cm
	Weight Measurement Range	0.2-30kg
	Height Measurement Tolerance	±5mm
	Weight Measurement Tolerance	0.2kg-5kg: ±5g; 5kg-20kg: ±10g; 20kg-30kg: ±15g
	Data Output Interface	RJ45*1; USB*2
	Data Transmit Methods	RJ45; WIFI; Bluetooth
	Bluetooth	5.2
	WIFI	WIFI 2.4
Others	Voltage	100V - 240V, 50/60Hz
	Rated Power	60VA
	Working & Storage Environments	Working Temperature 5°C – 40°C Working Humidity 10% – 95%RH (no condensation) Storage Temperature –10°C – 60°C Storage Humidity 10% – 95%RH (no condensation)
	Draduot Dimension	Atmospheric Pressure 70kPa – 106kPa
	Product Dimension	1279.61(L)*639.26(W)*397.5(H)mm
	Weight	about 24.2kg

Use Cases













YeoMed

- Shenzhen Yeo Medical Technology Co.,Ltd.
- Address: Floor 6, Building A, Wanhe Science and Technology Building, No.7 HuitongRoad, Guangming District, Shenzhen China.
- Factory Address: Room 402, Building C1, China Merchants Smart Park, No. 459 Qiaokai Road, Guangming District, Shenzhen China.
- Website: www.yeomedical.com